

CAD/CAM Orthotic Prescription and Design

Session 3

Total Contact Orthotic & Shoe Modifications for High Risk Feet

Presented By
Paul Graham

StepForce Training
12/09/2020

1

All about perspective



Copyright © Paul Graham 2020

2

DIABETICMed
DOI: 10.1111/j.1464-5491.2009.02835.x

Original Article: Complications

Plantar pressures in diabetic patients with foot ulcers which have remained healed

T. M. Owings, J. Apelqvist*, A. Sternström†, M. Beckert, S. A. Bus*, A. Kasperk, J. S. Ulbrecht† and P. R. Cavanagh**

Department of Biomedical Engineering, Center for Research in Assistive Devices, Cleveland Clinic, Cleveland, Ohio, USA; *Department of Endocrinology, University Hospital of Västerås, Västerås, Sweden; †Department of Orthopedic Surgery, University of Utah, Salt Lake City, Utah, USA; ‡Department of Endocrinology, University of Michigan, Ann Arbor, Michigan, USA; §Department of Orthopedics, University of Pennsylvania, Philadelphia, PA, USA; ¶Department of Orthopedics, Academic Medical Center, University of Amsterdam, Amsterdam, The Netherlands; **Foot Health, University of North Carolina at Chapel Hill, Chapel Hill, NC, USA; Accepted 1 August 2009

Abstract
Aims The recurrence of foot ulcers is a significant problem in people with diabetic neuropathy. The purpose of this study was to measure in-shoe plantar pressures and other characteristics in a group of neuropathic patients with diabetes who had prior foot ulcers which had remained healed.
Methods This was an epidemiological cohort study of patients from diabetes clinics of two Swedish hospitals. From a database of 2462 diabetic patients, 190 neuropathic patients with prior plantar ulcers of the foot (heel, ball or metatarsal heads) caused by repetitive trauma were identified and 49 patients agreed to participate. In-shoe and in-shoe plantar pressures were measured during walking. Data on foot deformity, activity profiles and self-reported behaviour were also collected.
Results Mean baseline plantar peak pressure at the prior ulcer site (164 kPa) was lower than in other published series, although the range was large (107–1152 kPa). Mean in-shoe peak pressure at the location averaged 207 kPa when measured with carbon-fibre sensors. Baseline peak pressure only predicted 3–5% of the variance of in-shoe peak pressure, indicating variation in the efficacy of the individual footwear prescriptions, primarily extra-depth shoes with custom insoles.
Conclusions We propose that the mean value for in-shoe pressures reported in these patients be used as a target in footwear prescriptions for patients with prior ulcers. Although plantar pressure is only one factor in a multifactorial strategy to prevent ulcer recurrence, the quantitative focus on pressure reduction in footwear is likely to have beneficial effects.
Diabet. Med. 26, 1141–1146 (2009)
Keywords baseline pressure measurements, in-shoe pressure measurement, neuropathy, threshold

Owings, T. M., et al. (2009). "Plantar pressures in diabetic patients with foot ulcers which have remained healed." *Diabet Med* 26(11): 1141-1146.


StepForce Training
12/09/2020
Copyright © Paul Graham 2020

3




van Netten, J. J., et al. (2018). "Diabetic Foot Australia guideline on footwear for people with diabetes." Journal of foot and ankle research 11(1): 2.


4



5



6




7



8




9



10




11



12




13



14



15



16



17



18



19



20



21

Footwear Modification

foot power®

what a relief!

Unit 11, 818 Pittwater Road (Enter via Delmar Parade), Dee Why NSW 2099
02 9972 4488

Step Force Training
12/09/2020 Copyright © Paul Graham 2020

22

Case Study

HISTORY

- 68-year-old women in overall good health
- Spinal Fusion - 2014; L4 & L5 – successful, but has resulted in numbness in the most of the foot up to the ankle. When barefoot she loses balance and struggles to walk. The Left foot has slight numbness through the forefoot.
- Right Foot - 2005. Fusion of the 2nd met and intermediate cuneiform - successful until 2 years ago. CT exam shows significant progression of osteoarthritis

Step Force Training
12/09/2020 Copyright © Paul Graham 2020

23

Case Study

HISTORY

- Metatarsal Fusion, left foot - 2019, 2,3 & 4 at midfoot. Also straightening and pinning L2/ - healing well but she is having repeated stress fractures due to the osteoporosis
- Has been advised that the R/foot will require further surgery but is concerned about having the stress fractures occur in the R/foot also.

Step Force Training
12/09/2020 Copyright © Paul Graham 2020

24

Case Study

2. EXAMINATION RESULTS

- Factors affecting the Force Pathway
 - Moderate High arch structure
 - Short first Metatarsal
- Compensation Available
 - Stiffness and pain in midfoot and L/hallux dorsiflexion
 - Lower leg Muscle strength shows weakness and imbalance
 - Balance and proprioception decreased ? Due to lumbar fusion
 - Tightness in the leg posterior muscle and soft tissues; L > R

StepForce
© 2019 Paul Graham 2020

25

Case Study



StepForce
© 2019 Paul Graham 2020

26

Orthoses don't work alone.






David Sutton CPedCM, AU
Clinical Advisor to Tekscan (Aust.) & Distributor.

27

The Why?




Van netten et al. 2018 paper produced for DFA. Guideline on footwear for people with Diabetes.


DFA is the peak national clinical and research body for diabetes-related foot disease in Australia.

12/09/2020

28





Van netten et al. 2018 paper produced for DFA. Guideline on footwear for people with Diabetes



Bus et al. Guidelines on the prevention of foot ulcers in persons with diabetes (IWGDF 2019 update). Diab Metab Res Rev. 2020. e3269

12/09/2020


29




“200kPa or reduce peak pressure by 30%...”

12/09/2020


30



31




32



33



34



35

Case presentation: Referred Patient, 74yr old male, 112 kg.
 Hallux rigidus RF, ulcer IP RF plantar. RL shorter 12mm, restricted ROM ankle.
 Ulcer is reoccurring over a 4 year period.
 Has Diabetes.
 Has retired from his volunteer job as a result.

MUSCLE STRENGTH TESTING (Oxford scale)


- KNEE FLEXORS: R 5 L 5
- KNEE EXTENSORS: R 4 L 5
- DORSIFLEXORS: R 4 L 4
- PLANTARFLEXORS: R 4 L 4
- EVERTERS: R 1 L 3
- INVERTERS: R 1 L 3

0/5 No contraction.
 1/5 Visible/palpable contraction, but no movement.
 2/5 Movement with no gravity resistance.
 3/5 Movement against gravity possible.
 4/5 Movement against gravity and some resistance.
 5/5 Movement against gravity with full resistance

36

We record 8 steps.

Sensor between the foot and the orthosis.




The image shows two foot pressure maps, one for the Left foot and one for the Right foot. A small sensor icon is visible on the medial ball of the foot in both. A vertical scale bar is positioned between the two maps.

37

First In-shoe pressure mapping (ISPM) assessment of the client.

Right Foot IP (ulcer site) 355kPa


Left Foot MPJ 344kPa



The image shows two foot pressure maps. The Left foot map has a peak pressure of 344 kPa at the MPJ. The Right foot map has a peak pressure of 355 kPa at the IP (ulcer site). A vertical scale bar is between the maps.

38

Peak pressure report




The report includes two line graphs and two foot pressure maps. The top graph is titled 'Peak Pressure versus Percent of Stance' and the bottom graph is 'Force versus Percent of Stance'. Both graphs show data for Left and Right feet. The foot maps show peak pressures of 344 kPa for the Left foot and 355 kPa for the Right foot.


39




40




41



42




43




44




45



46



47



48

Differential Table	3-Box - Foot Parameters - Difference Paired L versus R					
	Left pre intervention		L-R Difference		Right Post intervention	
COF Deviation (centimeters)	-0.6 to 0.9	0.3 to 0.9	-1.0 to -0.0	-0.7 to 0.6	0.3 to 1.5	-0.9 to -0.8
COF Excursion Index (%)	17%	7%	9%	15%	13%	1%
1st Peak (sec)	0.18	0.09	0.08	0.16	0.11	0.04
Trough (sec)	0.23	0.13	0.11	0.24	0.20	0.04
2nd Peak (sec)	0.44	0.30	0.14	0.47	0.45	0.02
Gait Curve 2-Peak Force Diff (kilograms)	31.29	34.43	-3.14	27.41	35.65	-8.24
Heel-Metatarsal Curves Crossing (sec)	-0.23	-0.27	0.04	-0.25	-0.23	-0.02
Heel Contact Time (sec)	0.64	0.71	-0.07	0.69	0.73	-0.04
Heel Maximum Force (%BW)	55%	30%	25%	53%	49%	4%
Heel Maximum Force (kilograms)	61.30	33.66	27.64	59.07	54.94	4.13

49

3-Box - Foot Parameters - Difference Paired L versus R						
Differential Table	Left pre intervention	Right	L-R Difference	Left Post intervention	Right	L-R Difference
COF Deviation (centimeters)	-0.6 to 0.9	0.3 to 0.9	-1.0 to -0.0	-0.7 to 0.6	0.3 to 1.5	-0.9 to -0.8
COF Excursion Index (%)	17%	7%	9%	15%	13%	1%
1st Peak (sec)	0.18	0.09	0.08	0.16	0.11	0.04

Normal Centre of Force excursion index for men (9.4 - 22.4) and for women (6.1 - 19.4)

Menz, H., Dufour, A., Riskowski, J., Hillstrom, H., & Hannan, M. (2013). Association of Planus Foot Posture and Pronated Foot Function With Foot Pain: The Framingham Foot Study. *Arthritis Care & Research*, 65(12), 1991-1999. doi: 10.1002/acr.22079

50

3-Box - Foot Parameters - Difference Paired L versus R						
Differential Table	Left pre intervention	Right	L-R Difference	Left Post intervention	Right	L-R Difference
COF Deviation (centimeters)	-0.6 to 0.9	0.3 to 0.9	-1.0 to -0.0	-0.7 to 0.6	0.3 to 1.5	-0.9 to -0.8
COF Excursion Index (%)	17%	7%	9%	15%	13%	1%
1st Peak (sec)	0.18	0.09	0.08	0.16	0.11	0.04
Trough (sec)	0.23	0.13	0.11	0.24	0.20	0.04
2nd Peak (sec)	0.44	0.30	0.14	0.47	0.45	0.02
Gait Curve 2-Peak Force Diff (kilograms)	31.29	34.43	-3.14	27.41	35.65	-8.24
Heel-Metatarsal Curves Crossing (sec)	-0.23	-0.27	0.04	-0.25	-0.23	-0.02
Heel Contact Time (sec)	0.64	0.71	-0.07	0.69	0.73	-0.04
Heel Maximum Force (%BW)	55%	30%	25%	53%	49%	4%
Heel Maximum Force (kilograms)	61.30	33.66	27.64	59.07	54.94	4.13

51

3-Box - Foot Parameters - Difference Paired L versus R


Differential Table	3-Box - Foot Parameters - Difference Paired L versus R					
	Left pre intervention	Right pre intervention	L-R Difference	Left Post intervention	Right Post intervention	L - R Difference
COF Deviation (centimeters)	-0.6 to 0.9	0.3 to 0.9	-1.0 to -0.0	-0.7 to 0.6	0.3 to 1.5	-0.9 to -0.8
COF Excursion Index (%)	17%	7%	9%	15%	13%	1%
1st Peak (sec)	0.18	0.09	0.08	0.16	0.11	0.04
Trough (sec)	0.23	0.13	0.11	0.24	0.20	0.04
2nd Peak (sec)	0.44	0.30	0.14	0.47	0.45	0.02
Gait Curve 2-Peak Force Diff (kilograms)	31.29	34.43	-3.14	27.41	35.65	-8.24
Heel-Metatarsal Curves Crossing (sec)	-0.23	-0.27	0.04	-0.25	-0.23	-0.02
Heel Contact Time (sec)	0.64	0.71	-0.07	0.69	0.73	-0.04
Heel Maximum Force (%BW)	55%	30%	25%	53%	49%	4%
Heel Maximum Force (kilograms)	61.30	33.66	27.64	59.07	54.94	4.13

52

3-Box - Foot Parameters - Difference Paired L versus R

Differential Table	3-Box - Foot Parameters - Difference Paired L versus R					
	Left pre intervention	Right pre intervention	L-R Difference	Left Post intervention	Right Post intervention	L - R Difference
COF Deviation (centimeters)	-0.6 to 0.9	0.3 to 0.9	-1.0 to -0.0	-0.7 to 0.6	0.3 to 1.5	-0.9 to -0.8
COF Excursion Index (%)	17%	7%	9%	15%	13%	1%
1st Peak (sec)	0.18	0.09	0.08	0.16	0.11	0.04
Trough (sec)	0.23	0.13	0.11	0.24	0.20	0.04
2nd Peak (sec)	0.44	0.30	0.14	0.47	0.45	0.02
Gait Curve 2-Peak Force Diff (kilograms)	31.29	34.43	-3.14	27.41	35.65	-8.24
Heel-Metatarsal Curves Crossing (sec)	-0.23	-0.27	0.04	-0.25	-0.23	-0.02
Heel Contact Time (sec)	0.64	0.71	-0.07	0.69	0.73	-0.04
Heel Maximum Force (%BW)	55%	30%	25%	53%	49%	4%
Heel Maximum Force (kilograms)	61.30	33.66	27.64	59.07	54.94	4.13

53




54



55



56



57
