

# Plantar Pressure Analysis & Diagnosis

USING PRESSURE TO PREDICT TREATMENT OUTCOMES

**StepForce Training**

Presented By Paul Graham

---

---

---

---

---

---

---

---

## Testing your treatment strategy

- Patient symptoms and history
- Asymmetries affecting the force pathway
- The body's ability to compensate
- Joint, muscle and soft tissue dysfunction
- Plantar Pressure data
  - dysfunction in pressure mapping, CoP and PvT graphs

**StepForce Training**

Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

## How interventions can really help



**StepForce Training**

Copyright © Paul Graham 2020

---

---

---


---

---

---

---

---




---

---

---


---

---

---

---

---




---

---

---

---

---

---

---

---




---

---

---

---

---

---

---

---

### Case Study 1. Plantar Heel Pain

#### 1. HISTORY

- 49-year-old mother of 3 in good overall health
- Present complaining of heel pain upon first weight bearing for last 2 months
- Has been working from home due to the virus and wearing flat slippers all day
- Noticed that the pain occurred around the time when she went out walking to improve her fitness and tried hill-walking to improve her cardio as she was 'feeling unfit' and was concerned about her weight gains that occurred from 'all the sitting'
- Her goals are to be pain free and able to walk regularly to lose weight



Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

### Case Study 1. Plantar Heel Pain

#### 2. EXAMINATION RESULTS

- No actual Injury and / or Surgical Factors
- Factors affecting the Force Pathway
  - Moderate Tibial varum
  - Moderate High arch structure
  - Short first Metatarsal
- Compensation Available
  - Some stiffness in midfoot and L/hallux dorsiflexion
  - Low leg Muscle strength is WNL for her age
  - Core Muscle strength is only fair and imbalanced
  - Tightness in the leg posterior muscle and soft tissues; L > R



Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

### Case Study 1. Plantar Heel Pain

#### 2. EXAMINATION RESULTS

- Adaptations
  - Dorsiflexed 1<sup>st</sup> metatarsals again L > R
  - Moderate FHL L/foot and Slight+ FHL R/foot
- Plantar Pressure Analysis
  - High loads sub L?/2<sup>nd</sup> and 3<sup>rd</sup> MTPs and Hallux IPJt
  - Central CoP with hesitations and minor blockages in forefoot
  - PvT graph shows L/arch elongation of 307Ms compared with 145 Ms in the R/foot



Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

### Case Study 1. Plantar Heel Pain

#### Plantar pressure of arch elongation




---

---

---

---

---

---

---

---

### Case Study 1. Plantar Heel Pain

#### Address Arch collapse



1. Use skin prep if possible
2. Apply 2.5 cm sports tape firmly
  - a) 1<sup>st</sup> piece to adduct forefoot
  - b) 2<sup>nd</sup> piece to plantarflex 1<sup>st</sup> Metatarsal
3. Unsure good adhesion
4. Start person walking over mat and record as soon as gait normalises

---

---

---

---

---


---

---

---

### Case Study 1. Plantar Heel Pain

#### Predicting Outcomes with Intervention




---

---

---

---

---

---

---

---

### Case Study 1. Plantar Heel Pain

#### Changing CoP trajectory



1. Tape the foot as to address arch collapse
2. Apply 5 cm sports tape firmly
  - a) 1<sup>st</sup> piece to evert forefoot
  - b) 2<sup>nd</sup> piece to support or Invert rearfoot
  - c) 3<sup>rd</sup> piece to midfoot if required
3. Unsure good adhesion
4. Start person walking over mat and record as soon as gait normalises

Step  
Force  
Training

Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

### Case Study 1. Plantar Heel Pain

#### Predicting Outcomes with Intervention



Step  
Force  
Training

Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

### Case Study 1. Plantar Heel Pain

#### Correct 1<sup>st</sup> MTP Joint function



1. Tape the foot as to address arch collapse
2. Add appropriate 5 cm tape to improve CoP trajectory; evert forefoot
3. Add 7mm felt pad under hallux, just distal to the 1<sup>st</sup> MTP joint
4. Unsure good adhesion
5. Start person walking over mat and record as soon as gait normalises

Step  
Force  
Training

Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

### Case Study 1. Plantar Heel Pain

#### Predicting Outcomes with Intervention



Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

### Case Study 1. Plantar Heel Pain

#### Predicting Outcomes with Intervention



Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

### Case Study 2. Lateral Forefoot Pain

#### 1. HISTORY

- 76-year-old man with a 50-year history of polio affecting his L/side
- L/leg shorter by 1.5 – 2cms with significantly weakened musculature
- Present complaining of severe pain in Left forefoot for some years
- Has tried physiotherapy, dry needling and other therapy that gives short term relief, but then "my calf cramps and my foot starts to hurt really hurting a lot".
- His goal is to have reduced pain and keep his mobility, which us scared of losing.

StepForce  
Training

Copyright © Paul Graham 2020

---

---

---

---

---

---


---

---



### Case Study 2. Lateral Forefoot Pain

#### Baseline CoP, PvT and Gaitline Analysis



Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

### Case Study 2. Lateral Forefoot Pain

#### Alter Pressure Versus Time Graphs



1. Decide how to alter the PvT Graph
2. Apply the pad under the heel that best suits the aim of the trial
3. Unsure good adhesion
4. Start person walking over mat and record as soon as gait normalises

StepForce  
Training

Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

### Case Study 2. Lateral Forefoot Pain

#### Baseline CoP, PvT and Gaitline Analysis



Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---



## Case Study 2. Lateral Forefoot Pain

### Addressing Overloading



1. Decide on the area that is overloaded
2. Decide how function can be improved to assist / resolve the overloading
3. Apply the pad under the heel that best suits the aim of the trial
4. Unsure good adhesion
5. Start person walking over mat and record as soon as gait normalises

Step  
Force  
Training

Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

## Case Study 2. Lateral Forefoot Pain

### Predicting Outcomes with Intervention



Step  
Force  
Training

Copyright © Paul Graham 2020

---

---

---

---

---


---

---

---

## Case Study 2. Lateral Forefoot Pain

### Predicting Outcomes with Intervention



Step  
Force  
Training

Baseline Assessment

With Interventions

Copyright © Paul Graham 2020

---

---

---

---

---

---

---

---

Questions?



Step  
Force  
Training

Copyright © Paul Graham 2020

---

---

---

---

---

---

---